Hybrid Knee R

Manufacture

Nabtesco Corporation

Accessibility Innovations Company
Assistive Products Department
35,Uozakihama-machi,Higashinada-ku
KOBE,658-0024,JAPAN
Phone : +81-78-413-2724
Fax : +81-78-413-2725
https://mobilityassist.nabtesco.com

Authorized representative for EU countries

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

















Nabtesco

Important Information

Cautions for handling Hybrid Knee R safely

Nabtesco Corporation (hereinafter referred to as “Nabtesco”) cannot foresee all potential re-sidual risks of Hybrid Knee R and risks resulting from human errors and the usage environment. Although there are many instructions and prohibitions for handling Hybrid Knee R (assem-bling, adjusting, and maintaining the prosthesis), all these matters cannot be described in this document or on the warning labels on the body of Hybrid Knee R.
Therefore, when handling Hybrid Knee R, it is necessary not only to observe the precautions stated in this document, but also to take the safety measures necessary for a prosthesis knee joint.
Particularly important matters concerning the safe handling of Hybrid Knee R are described below.
These matters apply to the persons who assemble and adjust Hybrid Knee R.
If you have a serious incident related to Hybrid Knee R, please report it to the manufacturer (contact information on the back cover) and the competent authority in your country.

Read this document thoroughly


Before handling Hybrid Knee R, thoroughly read this document, and sufficiently understand the contents.
In addition, please fully explain to users how to safely handle and care for this product, and be sure to follow the precautions in this manual.

1.Outlline of the Product

1-1 Function

● Specifications

Type	NK-3	NK-3SH
	Pyramid	Threaded head
Overall length	262.5mm	256.5mm
Dimension A	51mm	45mm
Dimension B	211.5mm	211.5mm
Dimension C	38mm	38mm
Weight	1,275g	1,285g
Knee flexion angle	140°	
Main Material	Titanium alloy, aluminum alloy	
Weight limit (specific conditions)	125kg (High Activity Level (K4), Hip Prosthesis: 100kg) Engineering Test Evaluation: ISO 10328 P6 (A 125kg)	
Activity level	K2 to K4 ※ Exclude activities at a competitive level	

* These specifications are subject to changes without prior notice.

● Lifetime: 5 years

● Basic configuration

The Hybrid Knee R is a knee joint that combines a hydraulic rotary damper and MRS system to control the stance phase, and a pneumatic cylinder to control the swing phase (see ❶).

EN

● Features

• Operational Principle

High stance stability and yielding function
The hydraulic rotary damper and unique MRS system reliably generate hydraulic resistance to prevent knee buckling. The yielding function also allows for safe descent of slopes and stairs.
When the prosthetic leg leaves the ground, the pneumatic cylinder activates, allowing the prosthetic leg to follow changes in walking speed. (See ❷)

• MRS System

When the knee joint is flexed, the hydraulic rotary damper produces hydraulic resistance, thus preventing abrupt buckling. The MRS (Mechanism of Reaction Force Sensing) system mechanically detects floor reaction force and sets hydraulic resistance to be ON or OFF. When floor reaction force exists at the heel side from the sensing point, hydraulic resistance is turned on. When it exists at the toe side, the hydraulic resis-tance is turned off. (See ❸)

• Yielding function

By having the heel land first, the patient can slowly flex the knee joint while applying body weight to the knee joint. After getting familiar with the yielding function, he/she can perform various acts such as descending a slope or stairs. (See ❹)

• Mode switching lever

Mode Switching Lever, you can switch the yielding mode, which is normally in normal mode, to safety mode. There are two modes, normal mode and safety mode, and you can switch between them by raising the lever to the top position. In safety mode, yielding is always active, preventing sudden knee flexion. It also allows you to sit down and descend slopes and stairs safely.
By switching to safety mode as needed, you can ensure safety and stability.
(See ❺（ **S safety mode, N normal mode**）)

1-2 Regarding use

● Intended purpose

Hybrid Knee R was designed and is manufactured for use as a prosthetic knee joint by above-knee amputated, knee disarticulated, and hip disarticulated patients. Do not use Hybrid Knee R for any other purposes. For hip prostheses, it is recommended to use a torsion adapter to prevent significant torsion from being applied to Hybrid Knee R. This device is intended for single patient multiple use.

● conditions of use

Do not use the Hybrid Knee R for any purpose other than as a prosthetic leg.
Assembly and adjustment must be performed by a prosthetist or orthotist.
Subcontracting this work to a third party is not permitted.

● Scope of application

It can be applied at a wide range of activity levels, including:

K2 Level : User can handle small environmental barriers such as curbs, steps, or uneven ground, both indoors and around the home.
K3 Level : User can handle most environmental barriers, and can walk at different speeds. In addition to simple walking, he/she can do light work and exercise as well.
K4 Level : User has physical abilities higher than basic walking. This includes children and athletes.
＊ Exclude activities at a competitive level

● Usable environmental conditions

Operating temperature range : -10℃ to +50℃
Storage temperature range : -20℃ to +60℃

EN

● Declaration of Conformity

Nabtesco Corporation hereby declares that the following Class I medical device complies with the essential health and safety requirements of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amend-ing Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC.

2. Safety Precautions

● **Definition of Symbols**

WARNING **Indicates a potentially hazardous situation which, if this symbol is not observed, could result in serious injury.**

CAUTION **Indicates a potentially hazardous situation which, if this symbol is not observed, could result in minor or moderate injury or property damage.**

● **Compatible Medical Devices (Prosthesis Components)**

Type	NK-3	NK-3SH
Proximal connectors	Female pyramidal connectors	M36 screw connectors
Distal connectors	φ 34 tube adapter	φ 34 tube adapter

● **Structural durability**

Hybrid Knee R was tested for 3 million walking cycles with a load of 125 kg which corresponds to the average walking distance for 3 years. We will not assume liability for ageing or damage of the product due to long-term usage.

ISO10328-P6-125kg*)

*) The body mass limit is not to be exceeded. For specific conditions and limitations of use, see the manufacturer's written instructions on intended use.

WARNING

■ Always explain the precautions for use to the user.
Improper use may result in falls or injury.

■ Upon detecting any abnormal noise, looseness, or drop in hydraulic resistance, discontinue use and contact your local sales representative/dealer.Continued use despite a detected abnor-mality may cause damage of parts, leading to falling.

■ Never use this product outside of its specifications or disassemble or modify it. Also, do not intentionally subject it to strong external impacts.
Doing so may result in personal injury or damage to the product or parts, leading to a fall.

■ Do not use for a person who weighs over 125 kg.
If used for a person who weighs over 125 kg, damage of parts may occur, leading to falling. This, however, does not preclude loading and unloading of baggage, etc. oc-curring in daily life.
* For K4 individuals and hip prosthesis the weight limit is 100 kg.

■ Do not use if you do not understand how to use the knee joint or the precautions for use.

■ Do not use if you have pain in the residual limb.

■ Do not use if you participate in sports that put a lot of strain on the knee joint.

■ When bending or straightening the knee joint, do not touch the moving parts. The hand can be caught in the swaying part of the pneumatic cylinder, causing serious injury.
*When standing up from a chair, never place the hand on the knee joint. Fingers can be caught between the knee joint plate and the hydraulic unit, causing serious inju-ry. If the patient needs assistance in standing up from a chair, we recommend placing hands on the armrests or on the seat.
*In addition to the above, placing a hand on the socket or placing both hands onthe sound leg can also help the patient stand up from a chair.
Provide instructions on the safe method according to the patient's circumstances.
(See ❻)

■ Never use this product for any purpose other than as a prosthetic knee joint.

EN

CAUTION

■ Adjustments should be made by a prosthetist using the adjustment method described in the instruction manual.
Incorrect adjustment can prevent pleasant walking.
Also, do not use the product without proper user training.
When training, be sure to ensure safety by holding onto the handrails, etc.

■ Do not allow contact with liquids such as water, salt water, chlorinated water, soapy water, gel soap, bodily fluids, and exudations. Do not use any detergent or solvent (thinner) for cleaning.
Doing so could result in rusting, discolor-ation, desiccation of the grease, resulting in a malfunction and abnormal noises.
Doing so may also result in poor appearance.

■ Use and store the product within the specified temperature range for environmental conditions. (See page 32.)
Otherwise, reduced functionality and deterioration or deformation of parts may occur, resulting in malfunction.

■ A socket should be formed and aligned so that neither the socket nor any other part contacts the pneumatic cylinder at the maximum flexion angle.
It will damage the pneumatic cylinder, im-peding normal walking. (See ❼)

■ Do not drop parts such as a screw in the frame.
Continued use after any part is dropped will damage the pneumatic cylinder, impeding normal walking.

■ Do not use beyond the expected service life.
Continuing to use the product when its functionality has deteriorated may result in a risk of falling.

■ If clothing or the foam cover gets caught in the moving parts, the clothing or foam cover may be torn.

■ Be careful in touching the thread head of the knee joint, which might have sharp edges and would injure your hand or fingers.

3. Before Use

Scope of Delivery

Please confirm that the following items are included.

- Hybrid Knee R 1
- Instruction Manual (this book) 1
- Extension Stopper Rubber 1

Assembly Procedure

Alignment

Perform static alignment according to the following procedure.

● Alignment of frontal plane

Complete alignment so that the load line passes through the center of the knee joint and falls in the heel center of the foot section.

● Alignment of sagittal plane

As shown in ❽ , complete alignment so that the weight load line passes through the range of 0 to 10 mm in front of the center of the knee axis of the knee joint. (See ❽)

CAUTION

When the device is not aligned as rec-ommended, excess load is imposed on the parts, accelerating wear. For a hip prosthesis, use of a torsion adapter is recommended.

Assembling the Tube (Pipe) (see ❾)

1. Oil and other substances can cause the pipe to rotate and slip, which can be dangerous. Be sure to degrease the pipe before installing it.
2. Loosen the tightening bolts on the outer clamp and insert the pipe from below, pushing it all the way in.
3. Align the split groove on the outer clamp with the split groove on the frame to install it.
4. Alternately tighten the two tightening bolts to a torque of 6.0 to 6.5 N m. Tightening one bolt will loosen the other, so alternately tighten them to the specified torque.

CAUTION

- **Always observe the torque values for tightening bolts. Tightening outside the specified torque values may result in damage or detachment of parts.**
- **For safe use, we recommend using pipes manufactured by Nabtesco. If you must use pipes from another manufacturer, make sure that the externaldimensional tolerance is a pipe diameter of Φ 34 + 0.05 , - 0.1mm.**
- **Do not use adhesive on tightening bolts.**

EN

4.How To Use

● How to use the yielding function


MRS characteristics

States changes depending on the position and timing of weight application, as shown in the diagram below.

(A) When body weight is applied to the heel, yielding function can work.

(B) When body weight is applied to the toe, yielding can not work.

(C) After sufficiently applying body weight to the heel, when the patient shifts body weight toward the toe while yielding function is retained, hydraulic resistance is sustained.

CAUTION

In order to make the yielding function effective,be sure to strike the floor with the heel first. Striking the floor with the toe first will turn off hydraulic resistance, meaning the body weight Stairs cannot be supported.
There is a risk that knee buckling may occur, causing falling-down. For safe use, fully understand the operational principle and proper walking method, and provide the patient with instructions.

● How to use the mode switching lever

Use the mode switching lever above the knee joint to select normal mode or safety mode.
(see ❺ (S , N))

1) Normal Mode (mode switching lever in lowest position) (See N)
This is the normal mode in which the MRS system switches yielding on and off based on the position of the ground reaction force.

2) Safety Mode (mode switching lever in highest position) (See S)
This mode ensures constant yielding.
This mode is useful when descending stairs or slopes, while seated, and when practicing yielding.

● Caution when using the mode switching lever

WARNING

- When switching the mode switching lever to safety mode, make sure it is raised all the way up and is in the highest position.
- To ensure safe use, please fully understand the yielding mechanism's operating principles and how to use it. Also, since the mode switching lever is operated by the wearer, please also instruct the wearer on how to operate it.

*The Hybrid Knee R is designed to assist the wearer's abilities, but it does not mean that anyone can easily descend steep slopes or stairs without training.

CAUTION

- Do not operate the Mode Switching Lever while loading with your toes.
 When the mode switching lever is operated during toe load, the yielding function does not activate, which can lead to an unexpected knee buckling.
 If you do so while loading with your toes, the lever may be difficult to turn, which could injure your fingers, or the lever may become deformed due to the load.
- Be careful when switching modes. Incorrect lever positions can cause knee buckling or stumbling, leading to falls. Ensure the lever is set properly before operating to avoid accidents.
- Do not subject the mode switching lever to impacts such as hitting it against nything or tipping it over from a leaning position with the knee joint. This could deform the mode switching lever and impair its switching function.
- If the mode swithcing lever gets caught on the foam cover or clothing and unexpected force is applied, it may rotate or be damaged.
- The rotation resistance of the mode switching lever may change with temperature changes.

EN

● About the foam cover

When wearing the foam cover, you can drill a hole in the area where the mode swiching lever would be and expose the mode swiching lever to reduce the risk of it getting caught on the foam cover or clothing.

● Precautions when handling prosthetic limbs

- The hydraulic rotary part will become hot with continued use of the prosthetic limb, so be careful not to get burned.
- Do not jump from high places or intentionally subject the prosthetic limb to strong impacts.
- Use the prosthetic limb under the specified environmental conditions.
- Do not allow the prosthetic limb to come into contact with liquids such as seawater, chlorinated water, sweat, or body fluids.
- Do not place the prosthetic limb in a high place or lean it unstable.
- Do not place objects on the prosthetic limb or step on it.
- Store the prosthetic limb under the specified environmental conditions.
 If using the prosthetic limb after leaving it unused for more than a month, bend and straighten the knee joint 2-3 times before use.


5.Adjustment

5-1 Adjusting the Stance Phase Control (hydraulic adjustment)

First, adjust the stance phase control. The hydraulic damper has two kinds of adjusting plugs, which are to be used as shown in Table 1. (see ⑩)

Table 1 How to Select Adjusting Plugs


	Adjusting plug	Effect	When to perform adjustment
Ⓐ	Resistance	Yielding strength (Hydraulic resistance)	Different body weight and activity level
Ⓑ	Sensitivity (OFF-timing)	Increasing effectiveness and sustainability (OFF-timing) of the yielding function.	When sticking occurs during swing transition, and when the yield duration is insufficient

The adjusting plugs have been set to the positions which are common to most people in the factory default setting.

If you can not find suitable positions, turn the plugs to their original position.

- **Ⓐ Hydraulic resistance adjusting plug: Set to the scale 10.**
- **Ⓑ Hydraulic sensitivity adjusting plug: Set to the scale 10.**

5-1-1 Adjusting the Hydraulic Resistance (see ⑪、⑫)

Insert a hex wrench (3mm) obliquely upward into the hydraulic resistance adjusting plug, and turn to change the hydraulic resistance.

- **Clockwise (Scale readings become smaller.) : Hydraulic resistance becomes higher at the stance phase.**
- **Counterclockwise (Scale readings become larger.) : Hydraulic resistance becomes lower at the stance phase.**

Adjust by using scale readings as a guide. Read the scale at the mark on the window.


EN

● Procedure

1. Make sure that the scale reads 10 and then have the patient stand within parallel bars.
2. Have the patient place his/her weight on the prosthesis and confirm that the yielding function works properly.
3. Have the patient feel two different levels of yielding strength by adjusting the scale as follows.
 - ① **Scale reading 0* (strongest) The reading "0" is where the plug is fully closed by turning clockwise.**
 - ② **Scale reading 15 (weak)**
4. Taking the reading 10 as the benchmark, finely adjust the adjusting plug and choose a desirable level of strength.

Guide for Intensity of Hydraulic Resistance

- Although hydraulic resistance varies by weight and activity, the standard level is generally between 5 and 15. If there are no problems with the strength of the hydraulic resistance, set it to 10.
- If your heel load is insufficient, the hydraulic resistance may feel weak.
 Please adjust your alignment, select the right foot, and receive walking instruction so that you can put enough weight on your heel.
 If your heel load is insufficient, yielding may not be sustained.

CAUTION


- **Do not use readings of 19 and above. When the plug is loosened one full turn or more from the fully closed position (scale 0) (past 19), a risk of oil leakage is created. Hydraulic resistance does not change even if the plug is further loosened.**

5-1-2 Adjusting Sensitivity (OFF-timing)

Adjust the sensitivity according to the needs shown in Table 1 (5-1 Adjusting stance phase control).
*Sensitivity adjustment is set to the standard setting when shipped from the factory. If there are no problems, no adjustment is necessary.
Use a hex wrench (3 mm) to turn the sensitivity adjustment plug on the right side as you face it.
(See ⓭、⓮)
When you turn it, the hydraulic resistance changes as follows:

- **Clockwise : The yielding function becomes more effective and more sustained.**
- **Counterclockwise : The yielding function becomes less effective and less sustained.**

The sensitivity adjusting plug can be adjusted up to one loosening turn away from the fully closed position. The factory default is set to the scale 10. If you can not find an appropriate position, it is recommended that you use the factory default setting.

CAUTION

- **Do not use readings of 19 and above.**
 The plug is designed for anti-dropping.
 Loosening the plug more than one turn from the fully closed position (scale 0) (past 19) may damage the parts.
- **Be careful not to tighten the MRS adjustment screw too much.**

EN

5-1-3 Walking on a Level Floor

1. **Have the patient walk within the parallel bars. Provide practice so that the patient fully extends the knee joint in the initial stance phase in normal walking on a level floor.**
2. **Confirm that there is no sticking caused by the yielding function at the time of transition from the stance phase to the swing phase.**
3. **Once the child becomes accustomed to walking, walk outside the parallel bars.**
 If there is catching, check the procedure and readjust it.

5-2 Swing phase control adjustment

To adjust swing phase control, use the two adjustment valves on the pneumatic cylinder to select the strength of flexion resistance and extension resistance. Adjust the extension adjustment valve and flexion adjustment valve on the pneumatic cylinder using the following procedure. First, adjust the flexion adjustment valve, and if terminal impact occurs during swing, adjust the extension adjustment valve. (See ⓯、⓰)

● Procedure

Flexion adjustment valve (⓱)

1. Use a 3mm hex wrench to turn the valve.
 - **Clockwise : Increases bending resistance**
 - **Counterclockwise : Decreases bending resistance**
2. Walk slowly while ensuring safety and make sure the knee joint bends in the optimal position. If it does not bend, turn it counterclockwise to decrease the resistance.
3. Gradually increase your walking speed and find the appropriate valve position.

Extension adjustment valve (⓲)

1. Use a 3mm hex wrench to turn the valve.
 - **Clockwise : Increases extension resistance**
 - **Counterclockwise : Decreases extension resistance**
2. Walk slowly while ensuring safety and make sure the knee joint is extended in the optimal position. If it is not extended, turn it counterclockwise to decrease the resistance.
3. Gradually increase your walking speed and find the optimal valve position.

If the extension resistance is too strong all at once, it can be dangerous as the knee will not be able to extend, so adjust it gradually.

Guidelines for swing phase control adjustment

- The adjustment range is one full turn counterclockwise from fully closed.
- The setting at the factory is one full turn counterclockwise from fully closed.
- Adjust to an appropriate position as this will vary depending on the strength of the wearer's swing.

CAUTION

- If the flexion resistance is set too high, there is a risk that the foot may unintentionally catch on the ground during the swing phase.
- If the extension resistance is set too high, the knee does not extend fully, thereby causing falling.
- Setting the extension resistance too high may affect the MRS system and prevent a smooth transition to the swing phase.
- To avoid abnormal air pressure resistance, adjust the flexion and extension valves gradually.
- Be careful not to tighten the valves too much.
 This may cause the valve to fall off, break, or make abnormal noises.

EN

5-3 Adjusting Yielding for Descending a Slope or Stairs

Before descending a slope or stairs, perform fine adjustments of hydraulic resistance and sensitivity.

Procedure

● **Adjustment for descending a slope**

- Adjustment of hydraulic resistance
 Select an appropriate strength of hydraulic resistance suitable for the angle of the declining slopes where the patient walks in his/her daily life.
- Adjustment of the sensitivity adjusting screw
 If sticking occurs at the toe-off when descending a slope or if the yielding function is not sustained, make fine adjustment of the sensitivity adjusting screw.

● **Adjustment for descending stairs**

- Adjustment of hydraulic resistance
 Select an appropriate strength of hydraulic resistance suitable for the stairs where the patient walks in his/her daily life.
- Fine adjustment of sensitivity
 Make fine adjustment of the sensitivity adjusting screw as needed.

How to Land on the Prosthesis When Descending Stairs

Instruct the patient to put the midfoot on the edge of the stairs. (See ⓳)

- Before using the prosthesis for walking on a slope or stairs, perform gait training. If appropriate use of yielding is deemed difficult after gait training, use of the yielding function should be prohibited. Failure to observe this instruction may cause a serious accident such as falling-down.

CAUTION

- If walking on a slope or stairs, the patient must be instructed to secure safety by holding a handrail. Due to the unfavorable floor conditions, there is a risk that the patient may lose balance, resulting in falling-down.
- The strength of the yielding effect varies depending on the wearer's weight, so please readjust it if there are significant changes in weight.

6. Maintenance parts

Replacing the Extension Stopper Rubber（See ⑳、㉑）

1. Bend the knee joint and remove the extension stopper rubber.

2. If there is any dirt or dust in the grooves of the knee joint, remove it. If there is any oil, degrease it.

3. Place the new extension stopper rubber with the protruding side facing up and the thicker side facing forward into the groove. Then press down firmly.

CAUTION **Do not use thinner or acetone to remove oil or grease.**

How to clean surface stains

If necessary, wipe off any surface dirt with a soft cloth dampened with fresh water.
*This product can be lightly washed with a small amount of fresh water, but be sure to immediately wipe off any moisture with a dry, soft cloth to ensure no moisture remains.

CAUTION **After wiping, be sure to remove any moisture or lint from the cloth.
Rust or entangled lint may cause abnormalities in the rotating parts, affecting bending and straightening of the knee joint.**

7. Disposal

The materials used in the main parts of the Hybrid Knee R are classified as follows:
When disposing of the product, please follow the rules of your local government.
Please instruct users not to dispose of the product as general waste at their own discretion.
Alternatively, if you return it to Nabtesco after notifying us, we will dispose of it for you.
(See ㉒)

	Part Name	Material classification	Note
1	Upper joint	metal	Titanium alloy
2	Knee Plate	metal	aluminum alloy
3	Extension stopper rubber	rubber	Nitrile rubber
4	Front Link	metal	aluminum alloy
5	Yielding switch lever	metal	stainless
6	Hydraulic Rotary	metal	Body: Aluminum alloy Rotor shaft: Iron alloy O-ring etc.: Rubber
7	Knee covers	plastic	Each material is specified
8	Fastening and supporting components	metal	iron alloy
9	Aluminum frame	metal	aluminum alloy
10	Air Cylinder	metal	Body: Aluminum alloy Rotor shaft: Iron alloy O-ring etc.: Rubber
11	Clamp	metal	stainless

8. Warranty Information

1) The warranty period for Hybrid Knee R should be confirmed with the Distributor. If failure or damage is caused by inadequacies of design or manufacture of the equipment supplied by us during the warranty period under normal usage, we will repair the relevant part(s).

2) The external part (extension stopper rubber) supplied with Hybrid Knee R is not covered by the warranty.

3) The warranty covers only defective items. We will not compensate for various expenses caused by the product failure or other damages.

4) In the following cases, the repair will be at the customer's expense even within the warranty period.

- ■ When failure or damage is caused by handling not conforming to the precautions stated in this document, dropping during use or giving impact not in the normal state of use.
- ■ When failure or damage is caused by incorrect use or improper modification or repair.
- ■ When failure or damage is caused by fire, earthquake, wind, flood, seawater, lightning or any other extraordinary natural phenomenon.
- ■ When failure or damage is caused by use of a device other than those designated by Nabtesco.
- ■ When failure or damage occurs after the warranty period has expired.

5) Liability

The manufacture will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document.

The manufacturer will not assume liability for damage caused by disregards of this document due to improper use or unauthorized modification, or for accident caused by the amputee's (user's) condition, especially due to the change of physical condition.

6) Resale and transfer of Hybrid Knee R

In thereby ensuring safe after-sales servicing, no reselling or assigning of the Hybrid Knee R is permitted.

7) Repair

- If the cause of defect is unclear, the measures will be determined through consultation between the customer and the Distributor.
- Repair of external painting will be performed at the user's expense. Ask the Distributor for estimation.
- If the knee joint must be repaired after the warranty period has expired, only replaced parts are warranted for a period of 6 manths from repairing.
- Hybird Knee R which has been used for 5 years after purchase may be unrepairable.

● Symbols Used on Labels

UDI Label (Packing Box)

Legal manufacturer

Made date

MD Medical device

REF Item number

GTIN Global Trade Item Number

SN Serial number

EC REP Authorized representative for EU countries


CE

Declaration of conformity according to the European Regulation 2017/745

Single patient multiple use

The maximum and minimum temperatures for storing, transporting, and using the product.

Body Mass Limit, device identification Label (Knee Joint Body)

TYPE : Product type

MAX. WEIGT : ~K3 : 125 kg K4 : 100 kg

You will find information about the K level on page 12.

MD SN CE Same as above

Tightening Torque Marking (Knee Joint Body) ❽ , ❾

6 ~ 6.5 Nm Tighten the clamp within this torque range.

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K2 Level : User can handle small environmental barriers such as curbs, steps, or uneven ground, both indoors and around the home.
K3 Level : User can handle most environmental barriers, and can walk at different speeds. In addition to simple walking, he/she can do light work and exercise as well.
K4 Level : User has physical abilities higher than basic walking. This includes children and athletes.
*Escludere attività agonistiche

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